

Hemochron® Response Coagulation Monitor



The Hemochron Response Coagulation Monitor is a whole blood coagulation monitor designed to be used at the point-of-care, it has an easy to read screen and has data management features

Features



Typical Manufacturer's Image

- Results displayed as whole blood, plasma equivalent, or INR for PT assay only
- Store 600 patient tests per well
- Store Operator and Patient Identification
- Allows supervisor to create and lock OID table
- Create and store up to 9 User Notes
- Attach up to 2 User Notes per record
- Download results to a PC
- Uses bar coded tubes
- Internal printer
- 2 level Liquid Quality Control lockout
- 2 level Electronic System Verification lockout
- Stores 300 QC results per well
- Test Menu
- Celite® Activated Clotting Time (ACT)
- Kaolin ACT
- Glass Activated ACT
- Activated Partial Thromboplastin Time (APTT)