



The ACT II automated coagulation timer is a portable point-of-care instrument used for hemostasis management that enables coagulation testing in cardiovascular and vascular surgery, catheterization laboratories, critical care units, and hemodialysis units. The ACT II provides multiple testing capabilities including various Activated Clotting Time (ACT) tests, the High Range Heparinase (HRHTC) test, the Prothrombin Time (PT) test and the General Purpose Cartridge (GPC) that can be used for the Activated Partial Thromboplastin Time (aPTT) test. Comes with DPU-414 Thermal Printer.

