

## GE HiSpeed Fx/i - Single Slice



Typical Manufacturer's Picture

HiSpeed Fx/i System includes:

- Large 70 cm gantry aperture
- Up to 350 mA Possible
- Plus or Minus (+/-) 30 degree tilt
- Scan time 1.0, 1.5, 2.0,3 sec
- 120 Sec. Spiral
- 1,2 or 3,5 & 10 mm Slice thickness
- Patient Positioning with Laser light
- Hi-resolution Hilite detector
- 0.7 sec partial scanning.

Optional Software features:

- Smart Prep
- Smart mA
- Smart Helical
- 3D Imaging Option Software
- MIP Option Software
- Real Time Reformation Software
- PAN / ZOOM option Software
- Display in Display Option Software
- Navigator



E-114 GIDC – Savli, Manjusar-Savli road,  
Dist: Vadodara, Tal: Savli, Pin Code 391775,  
Phone# +91-2667-26473, Fax# +91-2667-264736 Email: [info@somatechindia.com](mailto:info@somatechindia.com)  
Rev 10.10