

## GE HiSpeed Nx/I - Dual Slice



Typical Manufacturer's Picture

---

HiSpeed Nx/i System includes:

Large 70 cm gantry aperture  
Up to 350 mA Possible  
Plus or Minus (+/-) 30 degree tilt  
Scan time 0.8, 1.0, 1.5, 2.0,3 sec  
120 Sec. Spiral  
1,2 or 3,5 & 10 mm Slice thickness  
Patient Positioning with Laser light  
Hi-resolution Hilite detector  
0.7 sec partial scanning.

**Optional Software features:**

Smart Prep  
Smart mA  
Smart Helical



E-114 GIDC – Savli, Manjusar-Savli road,  
Dist: Vadodara, Tal: Savli, Pin Code 391775,  
Phone# +91-2667-264736, Fax# +91-2667-264736 Email: [info@somatechindia.com](mailto:info@somatechindia.com)

Rev 10.10

3D Imaging Option Software  
MIP Option Software  
Real Time Reformation Software  
PAN / ZOOM option Software  
Display in Display Option Software  
Navigator

Soma Tech Pvt. Ltd.



E-114 GIDC – Savli, Manjusar-Savli road,  
Dist: Vadodara, Tal: Savli, Pin Code 391775,  
Phone# +91-2667-264736, Fax# +91-2667-264736 Email: [info@somatechindia.com](mailto:info@somatechindia.com)  
Rev 10.10