


GE HiSpeed NX/i Pro – Dual Slice



Typical Manufacturer's Picture

The HiSpeed NX/i System with Pro Performance Package is a Premium Helical Whole Body Computed Tomography Scanner. The System Configuration is Designed to Produce Up to Two Images of Optimum Image Quality Per Rotation, Fully Simultaneous Scan Reconstruction, and Display.

- Dual-Slice Detector
- Interface: 1-Screen, Mouse Driven, Icon Based
- Tube heat capacity: 6.3 MHU
- Fastest rotation: 0.7 seconds
- Maximum scan technique: 440 mA at 120 kV
- kV settings: 80, 120, 140
- Nominal slice thicknesses (mm): 1, 2, 3, 4, 5, 6, 7, 10 (all 1X or 2X)



E-114 GIDC – Savli, Manjusar-Savli road,
Dist: Vadodara, Tal: Savli, Pin Code 391775,
Phone# +91-2667-264736, Fax# +91-2667-264736 Email: info@somatechindia.com
Rev 10.10